



Email : sales@laboratoryinstrumentindia.com

Product Name :
Polypropylene Shoe Covers

Product Code :
LBNY-00017-11000214



Description :

Anti-skid soles provide better traction to help avoid accidental slips
The breathable, lightweight fabric helps keep feet cool and comfortable as you work long hours.
Low-cost protection against dust, dirt, grime, and other dry particulates
These blue polypropylene shoe covers protect against dirt, grime, dust, and other particulates in non-hazardous environments. Anti-skid soles provide better traction to keep you on your feet and avoid accidental slips.
Available in standard and anti-skid options
Possible applications include food processing, pharmaceuticals, environmental cleanup, contract work, moving services, and more.

Technical Specification :

Polypropylene Shoe Covers



Laboratory Instrument India