



Email : sales@laboratoryinstrumentindia.com

Product Name :
Flat Bouffant Caps

Product Code :
LBNY-00017-11000215



Description :

Designed to protect against dirt, grime, and particulates in non-hazardous environments
Provide a particulate barrier and helps prevent hair from contaminating a product or work environment
Flat style
Latex-free elastic
Flat bouffant caps provide protection against dirt, grime, and particulates in non-hazardous environments.
Available in blue or white in a variety of sizes.
Low-linting, lightweight cover restrains hair and offers great ventilation
Possible applications include food processing, construction, pharmaceuticals, environmental cleanup, and more.

Technical Specification :

Flat Bouffant Caps



Laboratory Instrument India