



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Flat Bouffant Caps

**Product Code :**  
LBNY-00017-11000215



**Description :**

Designed to protect against dirt, grime, and particulates in non-hazardous environments  
Provide a particulate barrier and helps prevent hair from contaminating a product or work environment  
Flat style  
Latex-free elastic  
Flat bouffant caps provide protection against dirt, grime, and particulates in non-hazardous environments.  
Available in blue or white in a variety of sizes.  
Low-linting, lightweight cover restrains hair and offers great ventilation  
Possible applications include food processing, construction, pharmaceuticals, environmental cleanup, and more.

**Technical Specification :**

Flat Bouffant Caps



**Laboratory Instrument India**