



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Safety Bouffant Caps, 24 in.

**Product Code :**  
LBNY-00017-11000216



**Description :**

Bouffant caps keep hair in place and help ensure a clean and hygienic work area.  
Polypropylene, latex-free  
Bouffant caps  
Size is 24 in. and available in a variety of colors.  
24 in. diameter  
Available in variety of colors.

**Technical Specification :**

Safety Bouffant Caps, 24 in.



**Laboratory Instrument India**