



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Gamma Guard Sterile Cleanroom Shoe Covers

**Product Code :**  
LBNY-00017-11000218



**Description :**

Disposable sterile cleanroom shoe covers are made from high MVTR, innovative fabric. Each pair is individually packaged to help prevent cross-contamination.  
Sterile-Don packaging protects against contamination during donning  
Sterile shoe covers are made of comfortable fabric and include non-skid soles to help prevent accidental slips as you work. Shoe covers protect against liquids and particulates, so your shoes stay pristine. Each shoe cover is low-linting and sterilized to the highest Sterility Assurance Level.  
Helps keep wearers cooler vs many competitive fabrics  
Tear-resistant fabric ensures your safety in demanding environments  
The low-linting fabric helps reduce the release of unwanted particulates in your work area while tunnelized elastic also contains particle shed. Each pair of shoe covers is individually bagged to help prevent cross-contamination. Each bag includes a red indicator dot for proof of sterilization.  
Bound seams.

**Technical Specification :**

Gamma Guard Sterile Cleanroom Shoe Covers



**Laboratory Instrument India**

