



Email : sales@laboratoryinstrumentindia.com

Product Name :
Clear Polyethylene Apron

Product Code :
LBNY-00017-11000220



Description :

Ideal for food handling and preparation and general maintenance applications.
Lightweight embossed 20mil polyethylene apron with ties at neck
Easy to clean
Individually packaged
Economical protection for light- to medium-duty applications
Excellent resistance to chemicals, fats, oils and grease
Die cut for consistent sizing.

Technical Specification :

Clear Polyethylene Apron



Laboratory Instrument India