



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
High-Tech Conversions Non-Skid Cleanroom Shoe Covers

**Product Code :**  
LBNY-00017-11000221



**Description :**

Non-skid soles for enhanced traction  
Ideal for lab and light industrial use  
Non-skid, non-woven shoe covers are ideal for use in laboratories and manufacturing environments.  
Low-linting nonwoven fabric helps filter particulates for contamination control  
For use in laboratories and manufacturing environments.  
Elastic band gives secure yet comfortable fit around the shoes.

**Technical Specification :**

High-Tech Conversions Non-Skid Cleanroom Shoe Covers



**Laboratory Instrument India**