



Email : sales@laboratoryinstrumentindia.com

Product Name :
High-Tech Conversions Non-Skid Cleanroom Shoe Covers

Product Code :
LBNY-00017-11000221



Description :

Non-skid soles for enhanced traction
Ideal for lab and light industrial use
Non-skid, non-woven shoe covers are ideal for use in laboratories and manufacturing environments.
Low-linting nonwoven fabric helps filter particulates for contamination control
For use in laboratories and manufacturing environments.
Elastic band gives secure yet comfortable fit around the shoes.

Technical Specification :

High-Tech Conversions Non-Skid Cleanroom Shoe Covers



Laboratory Instrument India