



Email : sales@laboratoryinstrumentindia.com

Product Name :
Safety Pleated Polypropylene Bouffant Cap

Product Code :
LBNY-00017-11000224



Description :

Comfortable, pleated bouffant caps feature an encapsulated elastic headband to keep hair in place. Available in a variety of sizes and colors.

Disposable design for fresh, hygienic use every time.

Flexible, pleated design that stretches over the workers head

Polypropylene

100% Latex Free

Economical Protection

Sizes Available: 21, 24, 28 in

Technical Specification :

Safety Pleated Polypropylene Bouffant Cap



Laboratory Instrument India