



Email : sales@laboratoryinstrumentindia.com

Product Name :
Safety Isolation Gown

Product Code :
LBNY-00017-11000225



Description :

Durable, multi-use, low-cost isolation gowns provide economical protection.
Non-woven, spunbond polypropylene material is breathable, lightweight and low-lint.
Ideal for food manufacturing, light industrial, medical or anywhere a low-particulate barrier garment is needed.
Mildew- and flame-resistant
Latex-free and non-allergenic
One size fits all up to 3X-Large
Protect against dirt, grime and non-hazardous dry particulates
Rear entry and extra-long ties
Color: Yellow.

Technical Specification :

Safety Isolation Gown



Laboratory Instrument India