



Email : sales@laboratoryinstrumentindia.com

Product Name :
Safety SMS Lab Coat

Product Code :
LBNY-00017-11000226



Description :

Ideal for class 1000 environmental and above

Non linting

100% latex free and flame resistant

Individually inspected for stitching, sizing, material and workmanship

Commonly used in electronic, medical device and pharmaceutical industries.

Keystone SMS material is soft and durable, which makes it easy to work with. It offers excellence barrier protection while keeping the worker cool, comfortable and does not sacrifice mobility and comfort.

Technical Specification :

Safety SMS Lab Coat



Laboratory Instrument India