



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hi-Vis Vests

Product Code :
LBNY-00017-11000227



Description :

Breathable polyester mesh fabric

Hook and loop closure

Side elastic for improved comfort and fit

Ideal for additional levels of safety for low risks activities where ANSI compliance is not required, such as shopping cart retrievers, dog walkers, parking lot attendants, joggers, delivery vehicle drivers, recreational workers.

Available in hi-vis lime yellow and hi-vis orange.

Technical Specification :

Hi-Vis Vests



Laboratory Instrument India