



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Concealed Snap Front Cotton Lab Coats

**Product Code :**  
LBNY-00017-11000229



**Description :**

Flame-resistant coats ideal for manufacturing industry or in university labs  
Concealed snap front  
Side vent openings  
Flame-resistant  
Fabric composition: 88% Cotton/12% Nylon  
Lapel collar  
Hemmed sleeve ends  
One chest and two hip pockets.

**Technical Specification :**

Concealed Snap Front Cotton Lab Coats



**Laboratory Instrument India**