



Email : sales@laboratoryinstrumentindia.com

Product Name :
Concealed Snap Front Cotton Lab Coats

Product Code :
LBNY-00017-11000229



Description :

Flame-resistant coats ideal for manufacturing industry or in university labs
Concealed snap front
Side vent openings
Flame-resistant
Fabric composition: 88% Cotton/12% Nylon
Lapel collar
Hemmed sleeve ends
One chest and two hip pockets.

Technical Specification :

Concealed Snap Front Cotton Lab Coats



Laboratory Instrument India