




Email : sales@laboratoryinstrumentindia.com

Product Name :
Safety Cleanroom Compatible Polypropylene Lab Coats

Product Code :
LBNY-00017-11000230



Description :

Non-woven 100% Spunbond polypropylene
Strong, high tensile strength
Lightweight, breathable
Polypropylene lab coats are made of a breathable fabric that provides economical protection against dry particulates, dirt, and grime.
Cleanroom compatible
Low-cost solution to protect workers against dry particulates.

Technical Specification :

Safety Cleanroom Compatible Polypropylene Lab Coats



Laboratory Instrument India