




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Safety Cleanroom Compatible Polypropylene Lab Coats

**Product Code :**  
LBNY-00017-11000230



**Description :**

Non-woven 100% Spunbond polypropylene  
Strong, high tensile strength  
Lightweight, breathable  
Polypropylene lab coats are made of a breathable fabric that provides economical protection against dry particulates, dirt, and grime.  
Cleanroom compatible  
Low-cost solution to protect workers against dry particulates.

**Technical Specification :**

Safety Cleanroom Compatible Polypropylene Lab Coats



**Laboratory Instrument India**