



Email : sales@laboratoryinstrumentindia.com

Product Name :
Sterile Coveralls

Product Code :
LBNY-00017-11000233



Description :

Made of a comfortable fabric that protects against various liquids and particulates
Coveralls are designed with tunnelized elastic that contains particle shed
Serged seams eliminate raw edges that can create and release additional particulate contamination
Sterile coveralls offer enhanced safety with attached hood and boots along with a storm flap over the zipper.
Individually packaged and sterilized to a SAL of 10⁻⁶
Low-linting fabric helps reduce the release of unwanted particulates in the workspace
Storm flap over the zipper adds an extra layer of protection
Convenient thumb loops are included to keep sleeves from riding up
Ankle ties secure the attached boots
Attached hood and boots encourages quicker and easier donning processing
Elastic wrists and ankles help prevent particulates from entering or exiting the suit
Helps keep the wearer cooler.

Technical Specification :

Sterile Coveralls



Laboratory Instrument India

