



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Womens High-Performance Lab Coats

**Product Code :**  
LBNY-00017-11000234



**Description :**

One-piece, lined, notched lapel collar  
Left breast pocket, bar-tacked pencil stall, two lower pockets  
Side vent openings; back has princess seams  
With six-button closure  
80% polyester/20% combed cotton blend  
Pre-cure durable press finish  
White.

**Technical Specification :**

Womens High-Performance Lab Coats



**Laboratory Instrument India**