



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Polypropylene Bouffant Caps

**Product Code :**  
LBNY-00017-11000235



**Description :**

Bouffants caps  
Polypropylene  
Polypropylene bouffant caps are available in blue or white in a variety of sizes.  
Choose from blue or white color  
Packaged in quantities of 500 or 1000 caps per case.  
Available in 21, 24, or 28 in. sizes.

**Technical Specification :**

Polypropylene Bouffant Caps



**Laboratory Instrument India**