



Email : sales@laboratoryinstrumentindia.com

Product Name :
Protective Wear Unisex Cotton Lab Coat

Product Code :
LBNY-00017-11000236



Description :

Quick-release 5-snap front

Knit cuffs

Unisex Lab Coats are made of breathable, comfortable, 100% cotton fabric and feature a quick-release five-snap front and knit cuffs.

2 oversized waist pockets

Breast pocket with pencil slot

100% cotton.

Technical Specification :

Protective Wear Unisex Cotton Lab Coat



Laboratory Instrument India