



Email : sales@laboratoryinstrumentindia.com

Product Name :
Protective Wear Cleanroom ESD Scrub Pants

Product Code :
LBNY-00017-11000237



Description :

Controls static electricity
Elastic back and drawstring front for secure comfort
Cleanroom ESD scrub pants control static electricity while providing an ideal base-layer of protection for any cleanroom garment system.
Color-coded loops for easy size identification.
Elastic and snaps for leg opening.

Technical Specification :

Protective Wear Cleanroom ESD Scrub Pants



Laboratory Instrument India