



Email : sales@laboratoryinstrumentindia.com

Product Name :
Protective Wear Flame-Resistant Womens Lab Coat

Product Code :
LBNY-00017-11000238



Description :

Best in class flammability protection.
Safety sleeve cuff with tunneled elastic
Princess-belt on back
Reusable FR fabric
Quick-release snap front
Excellent heat-resistance properties.
FR identification patches on sleeves and back.
These flame-resistant lab coats were developed specifically to fit the female form while also providing the protection required in many of today's high-tech laboratories and research facilities.

Technical Specification :

Protective Wear Flame-Resistant Womens Lab Coat



Laboratory Instrument India