



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Aseptic Extra-Safe Lab Jackets — Light Pink

**Product Code :**  
LBNY-00017-11000240



**Description :**

Static-free — no clinging  
Made of soft three-layer SMS (Spunbond/Meltblown/Spunbond) fabrics.  
Autoclavable and washable up to five times  
Latex- free, fire retardant.  
Fluid- resistant, breathable  
Hip length  
Snap front, three pockets, knit cuffs and collar.

**Technical Specification :**

Aseptic Extra-Safe Lab Jackets — Light Pink



**Laboratory Instrument India**