



Email : sales@laboratoryinstrumentindia.com

Product Name :
Aseptic Extra-Safe Lab Jackets — Light Pink

Product Code :
LBNY-00017-11000240



Description :

Static-free — no clinging
Made of soft three-layer SMS (Spunbond/Meltblown/Spunbond) fabrics.
Autoclavable and washable up to five times
Latex- free, fire retardant.
Fluid- resistant, breathable
Hip length
Snap front, three pockets, knit cuffs and collar.

Technical Specification :

Aseptic Extra-Safe Lab Jackets — Light Pink



Laboratory Instrument India