



Email : sales@laboratoryinstrumentindia.com

Product Name :
Standard-Weight Polypropylene Shoe Covers

Product Code :
LBNY-00017-11000241



Description :

Elastic seam is turned in to reduce particle generation
With a non-marking, latex-based, anti-skid bottom surface
35g/m2 fabric is lightweight, durable, and breathable
Trimmed in elastic for a snug, adjustable fit
100% virgin spunbonded, non-linting polypropylene
Color: sky blue.

Technical Specification :

Standard-Weight Polypropylene Shoe Covers



Laboratory Instrument India