



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Standard-Weight Polypropylene Shoe Covers

**Product Code :**  
LBNY-00017-11000241



**Description :**

Elastic seam is turned in to reduce particle generation  
With a non-marking, latex-based, anti-skid bottom surface  
35g/m2 fabric is lightweight, durable, and breathable  
Trimmed in elastic for a snug, adjustable fit  
100% virgin spunbonded, non-linting polypropylene  
Color: sky blue.

**Technical Specification :**

Standard-Weight Polypropylene Shoe Covers



**Laboratory Instrument India**