



Email : sales@laboratoryinstrumentindia.com

Product Name :
PMD Pro safe Disposable Premium Economy Lab Coats,
Blue

Product Code :
LBNY-00017-11000242



Description :

Fluid resistant and breathable for protective comfort
Anti-Static treated with no clinging
Premium Quality durable unisex lab coats made of SMS offer flame retardant 3-layer SMS protection, fluid resistance, and are anti-static, and breathable.
Soft yet durable 3-Ply SMS latex-free material
Autoclavable or washable up to 5 times
Knit cuffs and collar for a snug feel
Snap closure front, and 3-pockets
Fire Retardant
Elongation testing showing 109 g/m² MD and 118.9 g/m² CD
Tested tensile strength of 119.7 g/m² MD and 61.6 g/m² CD
Knee Length.
Individually wrapped, sold in packs of 10.

Technical Specification :

PMD Pro safe Disposable Premium Economy Lab Coats, Blue

Laboratory Instrument India

