



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
PMD Pro safe Disposable Premium Economy Lab Coats,  
Blue

**Product Code :**  
LBNY-00017-11000242



**Description :**

Fluid resistant and breathable for protective comfort  
Anti-Static treated with no clinging  
Premium Quality durable unisex lab coats made of SMS offer flame retardant 3-layer SMS protection, fluid resistance, and are anti-static, and breathable.  
Soft yet durable 3-Ply SMS latex-free material  
Autoclavable or washable up to 5 times  
Knit cuffs and collar for a snug feel  
Snap closure front, and 3-pockets  
Fire Retardant  
Elongation testing showing 109 g/m<sup>2</sup> MD and 118.9 g/m<sup>2</sup> CD  
Tested tensile strength of 119.7 g/m<sup>2</sup> MD and 61.6 g/m<sup>2</sup> CD  
Knee Length.  
Individually wrapped, sold in packs of 10.

**Technical Specification :**

PMD Pro safe Disposable Premium Economy Lab Coats, Blue

---

# Laboratory Instrument India

