



Email : sales@laboratoryinstrumentindia.com

Product Name :
Shoe Cover with Conductive Strip

Product Code :
LBNY-00017-11000243



Description :

Latex based skid resistant bottom coating
Individually inspected for sizing, stitching and workmanship
100% virgin polypropylene into a nonwoven cloth-like single layer
Low lint cleanroom shoe cover. Shoe Cover with Conductive Strip features latex based skid resistant bottom.
Conductive Strips:
Constructed from a proprietary lint free co-polymer
Surface Resistance: $\geq 10^5$
Tension strength: 15 MPa.
Material Weight: 42gm/m²
Colors: Sky blue fabric, white bottom coating, black conductive strip.

Technical Specification :

Shoe Cover with Conductive Strip



Laboratory Instrument India

