



Email : sales@laboratoryinstrumentindia.com

Product Name :
Smooth Gowns Bulk Individually Bagged

Product Code :
LBNY-00017-11000244



Description :

Provide light-weight, impervious protection against bio-hazards. Smooth Gowns are in compliance with FDA's liquid and aerosol blood resistance testing requirements.
Easy, tear-away feature for quick gown removal and safe disposal
Contoured thumb loops secure gown sleeve for added protection
Redirects stress to the film material not the thermal seals to minimize breakage potential
Full-length open back affords better comfort and shoulder coverage
Provides liquid and aerosol blood resistance
Lightweight, non-linting construction.
Successfully tested for liquid and aerosol blood resistance under Synthetic Blood Penetration Resistance and Bacteriophage Penetration Resistance
In compliance with OSHA Bloodborne Pathogen Standard
Superior physical properties over imports
They have a light, non-linting construction for greater comfort.
In compliance with FDA requirements and tested in accordance with ASTM standards including Breaking Load, Tensile Strength and Elongation.

Technical Specification :

Smooth Gowns Bulk Individually Bagged

Laboratory Instrument India

