



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Smooth Gowns Bulk Individually Bagged

**Product Code :**  
LBNY-00017-11000244



**Description :**

Provide light-weight, impervious protection against bio-hazards. Smooth Gowns are in compliance with FDA's liquid and aerosol blood resistance testing requirements.  
Easy, tear-away feature for quick gown removal and safe disposal  
Contoured thumb loops secure gown sleeve for added protection  
Redirects stress to the film material not the thermal seals to minimize breakage potential  
Full-length open back affords better comfort and shoulder coverage  
Provides liquid and aerosol blood resistance  
Lightweight, non-linting construction.  
Successfully tested for liquid and aerosol blood resistance under Synthetic Blood Penetration Resistance and Bacteriophage Penetration Resistance  
In compliance with OSHA Bloodborne Pathogen Standard  
Superior physical properties over imports  
They have a light, non-linting construction for greater comfort.  
In compliance with FDA requirements and tested in accordance with ASTM standards including Breaking Load, Tensile Strength and Elongation.

**Technical Specification :**

Smooth Gowns Bulk Individually Bagged

---

# Laboratory Instrument India

