



Email : sales@laboratoryinstrumentindia.com

Product Name :

Aseptic Extra-Safe Lab Jackets — Ceil Blue

Product Code :

LBNY-00017-11000249



Description :

Autoclavable and washable up to five times
Latex- free, fire retardant.
Static-free — no clinging
Fluid- resistant, breathable
Snap front, three pockets, knit cuffs and collar
Made of soft three-layer SMS (Spunbond/Meltblown/Spunbond) fabrics.
Hip length.

Technical Specification :

Aseptic Extra-Safe Lab Jackets — Ceil Blue



Laboratory Instrument India