



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Aseptic Extra-Safe Lab Jackets — Ceil Blue

**Product Code :**  
LBNY-00017-11000249



**Description :**

Autoclavable and washable up to five times  
Latex- free, fire retardant.  
Static-free — no clinging  
Fluid- resistant, breathable  
Snap front, three pockets, knit cuffs and collar  
Made of soft three-layer SMS (Spunbond/Meltblown/Spunbond) fabrics.  
Hip length.

**Technical Specification :**

Aseptic Extra-Safe Lab Jackets — Ceil Blue



**Laboratory Instrument India**