



Email : sales@laboratoryinstrumentindia.com

Product Name :
Extra-Safe Lab Coats, Knee-Length — Cranberry

Product Code :
LBNY-00017-11000250



Description :

Autoclavable and washable up to five times
Latex-free, fire retardant.
Static-free — no clinging
Snap front, three pockets, knit cuffs and collar
Knee length
Made of soft three-layer SMS (Spunbond/Meltblown/Spunbond) fabrics.
Fluid-resistant, breathable.

Technical Specification :

Extra-Safe Lab Coats, Knee-Length — Cranberry



Laboratory Instrument India