



Email : sales@laboratoryinstrumentindia.com

Product Name :
Safety Black Reflective Safety Vest

Product Code :
LBNY-00017-11000251



Description :

Lightweight polyester solid material on front, mesh material on back
Black Reflective Safety Vest features include high visibility markings on the front and back, pockets, and a zipper front closure.

Clear front pocket

8 pockets: 5 outside, 3 inside

2 in. silver reflective material with 4 in. contrasting color trim

Badge/Microphone tabs located 3 in. below each shoulder seam.

Technical Specification :

Safety Black Reflective Safety Vest



Laboratory Instrument India