



Email : sales@laboratoryinstrumentindia.com

Product Name :
Meta Unisex Long Lab Coats

Product Code :
LBNY-00017-11000252



Description :

For men and women
Three patch-style pockets on outside, plus one inside pocket
Button-front closure
Two-piece swing belt
Made of 65% polyester, 35% cotton poplin
Length: 40 in. (102cm)
Washable
Not fluid resistant or antimicrobial. 4 pockets. Slits provide access to pants pockets.
Embroidery available upon request.
Tall and sizes up to 4X available upon request.

Technical Specification :

Meta Unisex Long Lab Coats



Laboratory Instrument India

