



Email : sales@laboratoryinstrumentindia.com

Product Name :

Protective Wear ESD Long Sleeve Scrub Top with Knit Cuff

Product Code :

LBNY-00017-11000253

**Description :**

Made with Sandwash material, this lightweight micro-denier fabric is woven from 100% continuous-filament polyester to eliminate any source for lint.

Widely used in pharmaceutical cleanroom operations

Pullover design with three-snap

ESD long sleeve scrub tops are made from non-linting, breathable fabric and are ideal base-layers in controlled environments.

Sleeves include color-coordinated filament knit cuff for extra comfort

5 in. side slits at garment bottom for added comfort and long-term durability

Henley-style collar that provides added room for easy on/off

Binding on collar helps maintain shape and uniformity

Raglan-style sleeves provide increased mobility

These scrubs provide an excellent base-layer of protection for any cleanroom garment system.

Technical Specification :

Protective Wear ESD Long Sleeve Scrub Top with Knit Cuff



Laboratory Instrument India

