



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Protective Wear ESD Long Sleeve Scrub Top with Knit Cuff

**Product Code :**  
LBNY-00017-11000253



**Description :**

Made with Sandwash material, this lightweight micro-denier fabric is woven from 100% continuous-filament polyester to eliminate any source for lint.  
Widely used in pharmaceutical cleanroom operations  
Pullover design with three-snap  
ESD long sleeve scrub tops are made from non-linting, breathable fabric and are ideal base-layers in controlled environments.  
Sleeves include color-coordinated filament knit cuff for extra comfort  
5 in. side slits at garment bottom for added comfort and long-term durability  
Henley-style collar that provides added room for easy on/off  
Binding on collar helps maintain shape and uniformity  
Raglan-style sleeves provide increased mobility  
These scrubs provide an excellent base-layer of protection for any cleanroom garment system.

**Technical Specification :**

Protective Wear ESD Long Sleeve Scrub Top with Knit Cuff



**Laboratory Instrument India**

