



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Shoe Covers

**Product Code :**  
LBNY-00017-11000254



**Description :**

Strong polypropylene and polyethylene material

Anti-skid

Waterproof

Fluid resistant and wear resistant

Shoe covers are made with a multi-layer fabric to provide water resistance and anti-skid properties and are ideal for labs and cleanroom environments.

Durable protection

Low linting, suitable for cleanrooms up to level ISO 5

Disposable

Available in Universal and XL sizes

Seamless sole and sides.

**Technical Specification :**

Shoe Covers



**Laboratory Instrument India**

