



Email : sales@laboratoryinstrumentindia.com

Product Name :
Poly Wear Smooth Disposable Gowns

Product Code :
LBNY-00017-11000255



Description :

Heat-sealed seam is behind sleeves—minimizes chance of contamination passing through.
Full-length wraparound with open back tailored to eliminate excess material
Tear-away collar for quick and easy removal
Lightweight, single-layer, low-density polyethylene gown for clinical and laboratory settings
Powder blue with smooth matte finish.
Thumb loop on sleeve keeps cuffs under gloves, enhancing leak protection
Gowns are folded for quick shake-out and easy donning.

Technical Specification :

Poly Wear Smooth Disposable Gowns



Laboratory Instrument India