



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Poly Wear Smooth Disposable Gowns

**Product Code :**  
LBNY-00017-11000255



**Description :**

Heat-sealed seam is behind sleeves—minimizes chance of contamination passing through.  
Full-length wraparound with open back tailored to eliminate excess material  
Tear-away collar for quick and easy removal  
Lightweight, single-layer, low-density polyethylene gown for clinical and laboratory settings  
Powder blue with smooth matte finish.  
Thumb loop on sleeve keeps cuffs under gloves, enhancing leak protection  
Gowns are folded for quick shake-out and easy donning.

**Technical Specification :**

Poly Wear Smooth Disposable Gowns



**Laboratory Instrument India**