



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Cryo-Apron

**Product Code :**  
LBNY-00017-11000257



**Description :**

Protection from splashes and contact cold in cryogenic and ultra-cold environments to -196°C.  
All aprons are 24 in. (61cm) wide  
Waterproof, lightweight  
Adjustable straps at neck and waist  
Color: blue.

**Technical Specification :**

Cryo-Apron



**Laboratory Instrument India**