



Email : sales@laboratoryinstrumentindia.com

Product Name :
Extra-Safe Lab Jackets — White

Product Code :
LBNY-00017-11000259



Description :

Static-free — no clinging
Autoclavable and washable up to five times
Latex-free, fire retardant.
Made of soft three-layer SMS (Spunbond/Meltblown/Spunbond) fabrics.
Snap front, three pockets, knit cuffs and collar
Hip length
Fluid-resistant, breathable.

Technical Specification :

Extra-Safe Lab Jackets — White



Laboratory Instrument India