



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Comfort Spunbonded Polypropylene Bouffant Caps

**Product Code :**  
LBNY-00017-11000260



**Description :**

Made from nonwoven fabric  
Latex-free  
Nonsterile  
Gentle synthetic elastic gives roomy, secure coverage without binding.  
Available in white or blue.

**Technical Specification :**

Comfort Spunbonded Polypropylene Bouffant Caps



**Laboratory Instrument India**