



Email : sales@laboratoryinstrumentindia.com

Product Name :
Class II High-Viz Vests

Product Code :
LBNY-00017-11000262



Description :

Two styles available
Surveyors Vest
Surveyor-style front pockets
Large arm holes
Zipper closure
100% polyester tricot
Two-tone striping with 2 in. wide (5.1cm) silver reflective tape
Color: Yellow
Standard Vest
Hook-and-loop closure
Machine-washable
Bright, lightweight polyester tricot fabric
2 in. wide (5.1cm) strip of reflective material—two 360° horizontal stripes and two vertical shoulder stripes
Color: Yellow.

Technical Specification :

Class II High-Viz Vests

Laboratory Instrument India

