



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Class II High-Viz Vests

**Product Code :**  
LBNY-00017-11000262



**Description :**

Two styles available  
Surveyors Vest  
Surveyor-style front pockets  
Large arm holes  
Zipper closure  
100% polyester tricot  
Two-tone striping with 2 in. wide (5.1cm) silver reflective tape  
Color: Yellow  
Standard Vest  
Hook-and-loop closure  
Machine-washable  
Bright, lightweight polyester tricot fabric  
2 in. wide (5.1cm) strip of reflective material—two 360° horizontal stripes and two vertical shoulder stripes  
Color: Yellow.

**Technical Specification :**

Class II High-Viz Vests

---

# Laboratory Instrument India

