



Email : sales@laboratoryinstrumentindia.com

Product Name :
High Performance SMS Lab Coats

Product Code :
LBNY-00017-11000263



Description :

Constructed of SMS fabric that is strong and durable, yet lightweight and resistant to tears and puncture
Features knit collar, left chest pocket and two side pockets with snap front and knit wrist
Low linting lab coats for use in critical and low particulate clean environments.
Impede the spread of viruses, germs, impurities and infections
Provides fluid and particulate safety barrier without sacrificing mobility and comfort.

Technical Specification :

High Performance SMS Lab Coats



Laboratory Instrument India