



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
High Performance SMS Lab Coats

**Product Code :**  
LBNY-00017-11000263



**Description :**

Constructed of SMS fabric that is strong and durable, yet lightweight and resistant to tears and puncture  
Features knit collar, left chest pocket and two side pockets with snap front and knit wrist  
Low linting lab coats for use in critical and low particulate clean environments.  
Impede the spread of viruses, germs, impurities and infections  
Provides fluid and particulate safety barrier without sacrificing mobility and comfort.

**Technical Specification :**

High Performance SMS Lab Coats



**Laboratory Instrument India**