



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Latex-Free Polypropylene Beard Cover

**Product Code :**  
LBNY-00017-11000264



**Description :**

Features a filament covered, unexposed elastic band for secure fit and closure  
Disposable SPP beard cover  
Manufactured with 100% non-linting, non-woven, spunbond polypropylene fabric  
Tested and manufactured in ISO compliant facilities under individual inspection and strict process  
Manufactured with 100% non-linting, non-woven, spunbond polypropylene fabric to provide comfort and protection. Recommended for Class 10,000 (ISO 7) cleanroom environments.  
Recommended for use in a cleanroom Class 10,000+ (ISO 7+) critical environment.

**Technical Specification :**

Latex-Free Polypropylene Beard Cover



**Laboratory Instrument India**