



Email : sales@laboratoryinstrumentindia.com

Product Name :
Latex-Free Polypropylene Beard Cover

Product Code :
LBNY-00017-11000264



Description :

Features a filament covered, unexposed elastic band for secure fit and closure
Disposable SPP beard cover
Manufactured with 100% non-linting, non-woven, spunbond polypropylene fabric
Tested and manufactured in ISO compliant facilities under individual inspection and strict process
Manufactured with 100% non-linting, non-woven, spunbond polypropylene fabric to provide comfort and protection. Recommended for Class 10,000 (ISO 7) cleanroom environments.
Recommended for use in a cleanroom Class 10,000+ (ISO 7+) critical environment.

Technical Specification :

Latex-Free Polypropylene Beard Cover



Laboratory Instrument India